WelchAllyn

9600 Plus Calibration Tester

© 2019 Welch Allyn. All rights are reserved. To support the intended use of the product described in this publication, the purchaser of the product is permitted to copy this publication, for internal distribution only, from the media provided by Welch Allyn. No other use, reproduction, or distribution of the publication, or any part of it, is permitted without written permission from Welch Allyn.

Welch Allyn assumes no responsibility for any injury to anyone, or for any illegal or improper use of the product, that may result from failure to use this product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual.

Welch Allyn is a registered trademark of Welch Allyn. Braun ThermoScan is a registered trademark of Braun Corporation.

Printed in Mexico

 729156

DIR 80020333 Ver B
Revised: 2019-07

Regulatory Affairs Representative
Welch Allyn Limited
Navan Business Park
Dublin Road
Navan, County Meath,
Republic of Ireland

Contents

Introduction	1
Symbols and Descriptions	2
Cautions and Warnings.	4
Unpacking the 9600 Plus Calibration Tester.	5
Controls and Indicators	7
LCD Display	7
Power On Sequence.	9
Front of the 9600 Plus Calibration Tester.	9
Operation.	13
9600 Plus Calibration Tester Setup.	13
Changing the Calibration Set Point Temperature	13
Welch Allyn SureTemp and SureTemp Plus Thermometers	15
Braun ThermoScan Pro 4000 Model.	17
Braun ThermoScan Pro 6000 model.	20
Troubleshooting inaccurate readings	25
Maintenance	27
Cleaning	27
Problem Solving	28
Specifications	31
Electromagnetic Compatibility	32
Warranty and Service	35
Limited Warranty	35

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1

Introduction

The 9600 Plus Calibration Tester verifies calibration of the following thermometers:

- Welch Allyn SureTemp[®] Plus models (690, 692 and OEM)
- Braun ThermoScan[®] Pro 4000 and Pro 6000

The 9600 Plus Calibration Tester provides a tightly regulated temperature source that eliminates the setup, precision calibration, and maintenance requirements of water baths and other general purpose calibration equipment. The user takes the “temperature” of the 9600 Plus Calibration Tester with the knowledge and assurance that precise readings are obtained because the unit is factory calibrated to National Institute of Standards and Technology (NIST) or equivalent devices. User adjustments are not required or recommended.

The “blackbody” of the 9600 Plus Calibration System acts as a heat reservoir for the thermistor probe and as a heat radiator for the infrared ear thermometers.

The unit contains solid-state circuitry for long-term stability with minimal maintenance. Welch Allyn uses only precision components in the construction of the 9600 Plus Calibration Tester. To maintain traceability to NIST or equivalent, return the 9600 Plus Calibration Tester to the factory for annual calibration.

Symbols and Descriptions

Double Insulation.

Direct Current.

Handle with Care.

Attention! Caution; Consult Directions for Use for more information.

Transport/Storage Temperature

Relative Humidity

Separate collection for electronic equipment

Recyclable

Manufacturer

SN

Serial Number

REF

Product Identifier

#

Reorder Number

EC REP

Authorized Representative in the European Community

Cautions and Warnings

A caution statement in this manual identifies a condition or practice which, if not corrected or discontinued immediately, could lead to equipment failure, equipment damage, or data loss.

Caution Do not immerse or soak the 9600 Plus Calibration Tester in any type of fluid.

Caution Only trained Biomedical personnel should operate the 9600 Plus Calibration Tester and conduct calibration verification of thermometers.

Caution When not using the Thermistor or Ear Device Ports, keep them covered with the Dust Cover.

Caution Store thermometers for testing in the same room as the 9600 Plus Calibration Tester for approximately 30 minutes prior to testing to allow the thermometer to reach ambient room temperature.

Caution Do not clean the Thermistor or Ear Device Ports. Doing so may damage the surface. Only a Welch Allyn authorized Service Center should clean the Device Ports.

WARNING

WARNING Use of this equipment is limited to the biomedical laboratory environment or equivalent.

WARNING Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING Portable RF equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the 9600 Plus Calibration Tester, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Unpacking the 9600 Plus Calibration Tester

Unpack the 9600 Plus Calibration Tester and inspect for missing items. Retain the shipping materials to return the device to Welch Allyn for recalibration.

The following items are enclosed:

- A 9600 Plus Calibration Tester
- A package containing a power supply and four mains plug inserts
- The Directions for Use
- A Carrying Case


If any of these items are missing after opening, please contact Welch Allyn Customer Service at:

<http://www.welchallyn.com/en/about-us/locations.html>

2 Controls and Indicators

LCD Display

Figure 1. LCD Display and Front Buttons

Calibration Set Point Temperature

There are three Calibration Set Point temperatures used to test the calibration of the thermometer. The three temperatures are 36.0° C (96.8° F), 38.5° C (101.3° F), and

41.0° C (105.8° F). The LCD displays the current Calibration Set Point Temperature in the upper left corner of the LCD.

Current Cavity Temperature

The current temperature of the 9600 Plus Calibration Tester is indicated by the large numbers in the center of the display (Figure 1).

As the device is transitioning to the current Calibration Set Point Temperature, the actual temperature flashes on the display. When the device reaches the Calibration Set Point Temperature, the temperature in the center of the display stops flashing.

Temperature Units

The 9600 Plus Calibration Tester is capable of displaying the temperature in degrees Fahrenheit or Celsius.

Set Point Mode

The LCD will display either a "D" (Default) or "L" (Last Used) in the upper left hand corner (Figure 1). The 9600 Plus Calibration Tester has two set point modes of operation. The Default mode (from the factory) causes the 9600 Plus Calibration Tester to heat up to the lowest set point (36.0° C / 96.8° F) upon startup. The Last Used mode causes the 9600 Plus Calibration Tester to heat up to the last temperature setting that was used.

To change the Set Point Mode, press and hold the C/F and Arrow buttons below the LCD display simultaneously until a single audio beep is heard. A "D" or an "L" will appear in the upper left corner of the display for the new mode.

Power On Sequence

The 9600 Plus Calibration Tester receives power from a wall-mount power pack that plugs into the back of the unit. There is no ON/OFF switch. On power up, the 9600 Plus Calibration Tester performs internal self-checks, displays an all-segments test pattern, beeps and displays the device's current temperature. The 9600 Plus Calibration Tester immediately applies heat, if needed, to reach the selected Calibration Set Point Temperature.

Front of the 9600 Plus Calibration Tester

Temperature Selection Button

This button allows the user to select one of three set point temperatures: 36.0°C (96.8°F), 38.5°C (101.3°F), or 41.0°C (105.8°F).

To scroll from one set point to the next, press and hold the Temperature Selection button until a beep is heard. The newly selected set point appears in the upper left corner of the LCD display. The device's current temperature is displayed, will start to flash, and continue flashing until the cavity reaches the equilibrium at the new set point.


Celsius/Fahrenheit Selection Button

This button allows the user to select degrees Celsius or degrees Fahrenheit. To change the temperature units, press and hold the C/F button until a beep is heard. The temperatures will appear in either "°C" or "°F", see [\(Figure 1\)](#).

Caution Verify the thermometer being tested and the 9600 Plus are set to use the same temperature scale before running a test. Reference the thermometers directions for use to change its temperature scale

Figure 2. Front of the 9600 Plus Calibration Tester

Dust Cover

This black, plastic cover keeps dust from entering the Thermistor and Ear Device Ports. To expose the Thermistor and Ear Device Ports, rotate the Dust Cover counter clockwise. To cover the ports, rotate the Dust Cover clockwise.

Caution When not using the Thermistor or Ear Device Ports, keep them covered with the Dust Cover.

Thermistor and Ear Device Ports

The Ear Device Port is the larger of the two ports located under the dust cover and is for the Braun ThermoScan Pro 3000 and Pro 4000 ear thermometers ([Figure 2](#)).

The Thermistor Device Port is the smaller port located under the dust cover and is for the SureTemp and SureTemp Plus thermometers.

The surface texture inside the Ear Device Port is painted with a special dull, black paint to enable Ear device readings. You can only see its surface under bright light. The Thermistor Device Port is not painted black. Keep both clean and dry.

The coating in the device ports of the unit is fragile, so please follow the cleaning instructions carefully. Do not touch the surface, and do not allow any foreign objects to be inserted into the ports. Avoid spilling liquids onto the unit. If liquid or other contaminants enter the unit, Welch Allyn can not assure proper calibration and you must return the unit to Welch Allyn for recalibration.


3 Operation

9600 Plus Calibration Tester Setup

Place the 9600 Plus Calibration Tester on a level surface away from sunlight, drafts, and other sources of heat or cold. The 9600 Plus Calibration Tester takes approximately 20 minutes to heat to the lowest setting. When testing several thermometers at all three temperatures, it is recommended to test all probes at one Calibration Set Point Temperature before proceeding to the next Calibration Set Point Temperature. To further expedite testing start at the lowest Calibration Set Point Temperature. The 9600 Plus Calibration Tester does not have an internal fan, causing the cool down time to be greater than the time it takes to heat up.

Changing the Calibration Set Point Temperature

To scroll from one set point to the next, press and hold the Temperature Selection button until a beep is heard. The newly selected set point appears in the upper left corner of the LCD display. The device's current temperature is displayed, starts to flash, and continues flashing until the cavity reaches the equilibrium at the new set point.

Figure 3. Power Adapter and Mains Plug Inserts

For the setup:

1. Choose the proper mains plug insert and slide it over the two prongs in the power converter (Figure 3).
2. Plug the power adapter into the 9600 Plus Calibration Tester (Figure 2) and the opposite end into a wall outlet.
3. Observe the Set Point Mode in the upper left hand corner of the LCD display (Figure 1). If the unit displays a "D", it is in Default Mode and will heat to the lowest Set Point Temperature. If you do not want to conduct testing at this Set Point Temperature, press and hold the Temperature Selection button to select the desired setting. The temperature display will flash before staying on continuously to indicate the 9600 Plus Calibration Tester has stabilized and is ready for use.


Welch Allyn SureTemp and SureTemp Plus Thermometers

Caution Store thermometers for testing in the same room as the 9600 Plus Calibration Tester for approximately 30 minutes prior to testing to allow the thermometer to reach ambient room temperature.

Note Prior to testing visually inspect the probe for any physical damage. Damaged probes must be replaced.

Figure 4. 9600 Plus Calibration Tester with the Welch Allyn SureTemp Plus Thermometer

To begin calibration verification of the thermometers:

1. Remove the probe from the probe well and clean it with either a 70% isopropyl alcohol solution, a 10% chlorine bleach

solution, or a non-staining disinfectant. Let the probe air dry. Do not apply a probe cover.

2. Place the thermometer in Monitor Mode, refer to the thermometer's Operator's Manual.
3. Insert the probe into the Thermistor Device Port ([Figure 4](#)).
4. Wait for approximately one minute or until temperature on the thermometer is stable for ten seconds. Compare the thermometer's temperature reading to the 9600 Plus' Calibration Set Point Temperature. If the temperatures are within $\pm 0.1^{\circ}\text{C}$ ($\pm 0.2^{\circ}\text{F}$), the thermometer is within calibration.
5. Test all available thermometers for calibration verification at the current Calibration Set Point Temperature. Proceed to the next Calibration Set Point Temperature, see "[Changing the Calibration Set Point Temperature](#)" ([page 13](#)).

Braun ThermoScan Pro 4000 Model


Caution Store thermometers for testing in the same room as the 9600 Plus Calibration Tester for approximately 30 minutes prior to testing to allow the thermometer to reach ambient room temperature.

Caution The ambient temperature must be stable and within the range of 18.3°C (65.0°F) – 26.7°C (80.0°F).

Caution Prolonged body contact with the thermometer may impact measurements. When instructed to wait before performing the next step place the thermometer on the work surface.

Figure 5. 9600 Plus Calibration Tester with the Braun ThermoScan Pro 4000

To begin calibration verification of the thermometers:

1. Check the probe tip lens window for displacement. There shouldn't be any gap in the seam formed between bezel and the lens window. If a gap exists the unit is damaged and must be replaced.
2. Gently wipe the probe tip with a cotton swab slightly moistened with isopropyl alcohol using care not to displace the lens window; remove excess alcohol with a clean cotton swab, and let air dry for 5 minutes. Do not use any chemical other than alcohol to clean the probe window.
3. Place the Braun ThermoScan Pro 4000 in Calibration Mode using the following steps.
 - a. Push and release the I/O MEM button to turn the product on. The display will show symbols and functions. The Pro 4000 will perform an automatic self check. After a few seconds you will hear a beep and see three dashed lines, a sound symbol, and ° C or ° F on the display.
 - b. Push the I/O MEM button again and keep it pressed for the following steps:
 - After approximately 3 seconds the "OFF" symbol will flash on the display (keep the button depressed).
 - When you hear a beep, release the button immediately.
 - The Pro 4000 is now in calibration check mode and the display is flashing and showing the "CAL" symbol.
4. After entering calibration mode (CAL), wait 1 minute before taking the first temperature.

Caution Prolonged body contact with the thermometer may impact measurements. When instructed to wait before performing the next step place the thermometer on the work surface.

5. Apply a new probe cover.
6. Place the probe firmly into the Ear Device Port (Figure 5).
7. Wait approximately 3 seconds, press the thermometer **Start** button, and watch for the "Exac Temp" light to flash.
8. Leave the thermometer in the 9600 Plus Calibration Tester until you hear a beep.
9. Remove the Pro 4000 Thermometer from the 9600 Plus Calibration Tester and read the temperature in the thermometer's display. If the temperatures are within $\pm 0.2^{\circ}\text{C}$ ($\pm 0.4^{\circ}\text{F}$), the thermometer is within calibration.
10. Wait 1 full minute before taking another reading with the same thermometer. Repeated measurements in short sequence may cause higher readings.

Caution Prolonged body contact with the thermometer may impact measurements. When instructed to wait before performing the next step place the thermometer on the work surface.

11. The probe cover must be replaced with a new probe cover before each reading, to remove any residual heat absorbed from the 9600 Plus Calibration Tester
12. The device will exit CAL mode after four minutes.
13. When using only one 9600 Plus Calibration Tester, test all available thermometers for calibration verification at the current Calibration Set Point Temperature. Proceed to the next Calibration Set Point Temperature, see "[Changing the Calibration Set Point Temperature](#)" (page 13).

Braun ThermoScan Pro 6000 model

Caution Store thermometers for testing in the same room as the 9600 Plus Calibration Tester for approximately 30 minutes prior to testing to allow the thermometer to reach ambient room temperature.


Caution The ambient temperature must be stable and within the range of 18.3°C (65.0°F) – 26.7°C (80.0°F)

Caution Set the temperature scales on the 9600 Plus to match the Braun ThermoScan PRO 6000's default temperature scales.

Caution Prolonged body contact with the thermometer may impact measurements. When instructed to wait before performing the next step place the thermometer on the work surface.

Caution On older versions of the 9600 Plus, the dust cover screw can interfere with the proper positioning of the Braun 6000 in the ear device port. Use the Braun 6000 test fixture (Welch Allyn part number 411690) to ensure proper positioning.

Figure 6. 9600 Plus Calibration Tester with the Braun ThermoScan Pro 6000

A Braun Pro 6000 Calibration Fixture (material number 411690) is available for use with older 9600 Plus calibration testers where the dust cover screw might interfere with the proper alignment of the Braun 6000 in the ear device port.

The older models have a dust cover screw that rises above the O-ring washer. Newer models have a dust cover screw that is recessed into the washer.

When using the calibration fixture, follow the instructions provided with the fixture to install it to the calibration tester.

To begin calibration verification of the thermometers:

1. Check the probe tip lens window for displacement. There should not be any gap in the seam formed between the bezel and the lens window. If there is a gap, the unit is damaged and must be replaced.

2. Gently wipe the probe tip with a cotton swab slightly moistened with isopropyl or ethyl alcohol. Use care not to displace the lens window. Remove excess alcohol with a clean cotton swab, and let air dry for 5 minutes. Do not use any chemical other than alcohol to clean the probe window.
3. Place the Braun ThermoScan Pro 6000 in Calibration Mode using the following steps.
 - a. With the thermometer in sleep mode, press and release **C/F**.
 - b. IMMEDIATELY (while all segments of the LCD are illuminated during self-test) press and hold **C/F** and **Mem**.
 - c. After 5 seconds the thermometer will emit a long beep. During the beep release **Mem** and **C/F**. The thermometer is now in Calibration Check Mode, indicated by **CAL** displayed on the screen for 0.25 seconds every second
4. After entering calibration mode (CAL), wait 1 minute before taking the first temperature.

Note Do not hold the thermometer while waiting. This allows the thermometer to remain at ambient room temperature.

5. Apply a new probe cover. If applied properly a short beep sounds and the flashing probe cover icon will disappear from the display.
6. Wait for the screen to flash three dashes (- - -) alternating with **CAL**.
7. If you are using the Braun 6000 test fixture, place the Braun Pro 6000 in the fixture as shown in the assembly instructions (material number 721533). If the Braun Pro 6000 adapter is

not used, place the probe firmly into the ear device port (Figure 6). It is important to ensure that the device is perpendicular in the port, directly aligned with the calibration tester sensor and not at an angle or tilted. Firmly push the device down into the well to insure perpendicularity and that the probe is fully seated in device port. Visually check to ensure proper alignment. Do not apply too much pressure once the device is firmly seated, as this may cause your hand to waiver and tilt or move the device during testing.

8. Wait 5 seconds, press the thermometer **Measure** button, and watch for the green ExacTemp light to flash.
9. Leave the thermometer in the 9600 Plus Calibration Tester until the ExacTemp light stops flashing and you hear a beep.
10. Remove the Pro 6000 Thermometer from the 9600 Plus Calibration Tester and read the temperature in the thermometer's display. If the temperatures are within ± 0.2 °C (± 0.4 °F), the thermometer is within calibration.

Note If temperatures are outside ± 0.2 °C, verify the test procedure was strictly observed and repeat the calibration check two more times at this temperature. If any two subsequent test are within the acceptable limits the test may be considered a pass. If all three readings are out of calibration, call Welch Allyn for service.

11. Replace the probe cover to reset the thermometer before taking the next measurement. Wait for the thermometer screen to flash three dashes (- - -) alternating with **CAL**

Note Do not hold the thermometer while waiting. This allows the thermometer to remain at ambient room temperature.

Note The probe cover must be replaced with a new probe cover before each reading, to remove any residual heat absorbed from the 9600 Plus Calibration Tester.

12. Wait one full minute before taking another reading with the same thermometer. Repeated measurements in short sequence may cause higher readings.

Note Do not hold the thermometer while waiting. This allows the thermometer to remain at ambient room temperature.

13. When using only one 9600 Plus Calibration Tester, test all available thermometers for calibration verification at the current calibration set point temperature before proceeding to the next calibration set point temperature. See "[Changing the Calibration Set Point Temperature](#)" (page 13).

14. The device will exit CAL mode after 10 minutes of inactivity, or can be forced to exit by holding both **C/F** and **Mem** for 5 seconds. After exiting Calibration Check Mode the thermometer will return to Sleep Mode.

Troubleshooting inaccurate readings

1. Make sure that the sensor window is clean and free of debris
2. Always use new probe covers when calibration testing Braun thermometers.
3. Make sure that the ear and non-contact thermometers are in calibration mode.
4. Make sure that the thermistor thermometers are in monitor mode.
5. If the measurement takes longer than 5 seconds, or an **ERR** message appears in the display, verify the thermometer is in calibration mode.
6. Make sure that the thermometer is not affected by any radiant or conductive heating or cooling.
7. If no other cause is found, allow the thermometer to sit untouched for 5 minutes before re-testing.

4 Maintenance

Cleaning

To clean the 9600 Plus Calibration Tester:

1. Unplug the electrical cord from the back of the unit.
2. Using a cloth slightly dampened with a solution of warm water and mild soap, clean the outside of the unit. Take extra care so that soap and water do not enter the two device ports.
3. Rinse the cloth with water and thoroughly wring it out.
4. Wipe off the soap and water from the 9600 Plus.
5. Let the 9600 Plus Calibration Tester dry for 15 minutes.
6. Plug the electrical cord into the back of the unit.

Caution Do not clean the Thermistor or Ear Device Ports. Doing so may damage the surface. Only a Welch Allyn authorized Service Center should clean the Device Ports.

Problem Solving

Error	Problem	Solution
E 	Calibration Tester cannot reach Set Point Temperature because room or device is too warm.	Check ambient temperature is within range. Allow 30 minutes for thermal accommodation prior to testing.
E 	Calibration Tester cannot reach Set Point Temperature because room or device is too cold.	Check ambient temperature is within range. Allow 30 minutes for thermal accommodation prior to testing.
Voltage	Power adapter voltage problem	Call Welch Allyn for replacement
E0.1	Broken Heater, not heating at all	Return to Welch Allyn for service
E0.2	Heater run away (heater does not turn off at appropriate time)	Return to Welch Allyn for service
E0.3	Thermistor circuit fail	Return to Welch Allyn for service
E0.5	RCAL resistor circuit fail	Return to Welch Allyn for service
E0.7	PTB circuit fail	Return to Welch Allyn for service
E0.9	Catastrophic internal failure	Return to Welch Allyn for service
E1.1	Device not calibrated	Return to Welch Allyn for service

9600 Plus Calibration Tester Checks

If several thermometers do not pass the calibration verification tests, check the 9600 Plus Calibration Tester to determine if it is out of calibration. Before sending the unit to an authorized Welch Allyn service center, complete the checks listed below.

1. Check the calibration date on the side of the 9600 Plus Calibration Tester. If the 9600 Plus Calibration Tester was calibrated over 12 months ago, do not use it. Send it to an

authorized Welch Allyn Service center immediately for calibration.

2. Ensure no one has tampered with the 9600 Plus Calibration Tester. If someone has tampered with the unit or the calibration label(s), the calibration settings and warranty are void. Send it to an authorized Welch Allyn Service center immediately for calibration.
3. Ensure the 9600 Plus Calibration Tester is not exposed to extreme temperatures that can damage the precision thermistor. Exposure to high heat can cause the resistance versus temperature characteristic of the unit to shift, causing the 9600 Plus Calibration Tester to become inaccurate. Send it to an authorized Welch Allyn Service center immediately for calibration.

Environmental Checks

1. Verify that there is no excessive moisture or steam in the air. High humidity (above 90%) can cause electrical leakages that can temporarily shift the device out of calibration. If the thermometers and/or the 9600 Plus Calibration Tester are exposed to such conditions, allow each device to sit at normal room temperature and humidity for at least 30 minutes before conducting a calibration test.
2. Avoid exposure to extreme temperatures. Thermal gradients in the thermometer or the 9600 Plus Calibration Tester can cause unpredictable effects. Verify the 9600 Plus Calibration Tester is within the specified ambient temperature range of 18.3° C (65.0° F) to 26.7° C (80.0° F). Verify the ambient room temperature is stable and there are no drafts around the 9600 Plus Calibration Tester.

Specifications

Power Requirements	Voltage: 100 - 240 V 50/60 Hz
Current	Heating: 480 mA After Warm-Up: 185 mA
Power Consumption	Heating: 2.4 W After Warm-Up: 0.92 W
Calibration Tester	Length: 254 mm (10.00") Width: 105 mm (4.13") Height: 85 mm (3.34") Weight (without Power Pack): 0.75 kg (1.65 lbs)
Power Pack (approximate dimensions)	Length: 86 mm (3.38") Width: 51 mm (2.00") Height (with plug inserted): 62 mm (2.44") Weight: 0.18 kg (0.40 lbs)
Environment	Storage Temperature: 0°C – 50°C (32°F – 122°F) Ambient Operating Temperature: 18.3° C to 26.7°C (65° F to 80° F) Operating Humidity: 20% to 95% RH. Non-Condensing
Operating Specifications	Controlled Temperature of 9600 Plus Low Setting: 36°C (96.8°F) Intermediate Setting: 38.5°C (101.3°F) High Setting: 41°C (105.8°F) Temperature Accuracy: ±0.05° C (0.1°F) Maximum Pressure/Altitude: ≤ 2000 meters.
Heating and Cooling Times Tested at an ambient of 23.9°C (75°F)	Room Temperature to Low Setting: Approximately 20 min. Low Setting to Intermediate Setting: Approximately 6 min. Intermediate Setting to High Setting: Approximately 6 min. Cool from High Setting to Low Setting: Approximately 25 min.
Installation Category II Pollution Degree II	

Electromagnetic Compatibility

This product has been tested and found to comply with all applicable requirements of council directive 2014/30/EU and its amendments. This product has been tested to EN 61326-1:2013 and found to comply with the performance criteria and requirements of the standard.

This product is intended for use in a Class B electromagnetic environment such as a biomed laboratory for equipment testing, maintenance, etc.

This product is not intended for use in the medical / patient environment. During testing, the product displayed consistent temperature during testing and did not deviate more than +/- 0.05 deg. C.

If operated under conditions outside the testing parameters:

- If the product experiences a shut down, it is expected that it can be restored without any errors;
- If the temperature deviates more than +/- 0.05 deg. C during an EMC event, the product is expected to auto-recover when the EMC interference stops.

This product can radiate radio frequency energy and if not used in accordance with the instructions, which may cause harmful interference to other devices in close proximity. However, there is no guarantee that interference will not occur in a particular setting. If this product does cause harmful interference with other devices, which can be determined by turning the product on /off, the user is encouraged to try to correct the interference by one of more if the following measures:

- Reorient or relocate the receiving device
- Increase the separation between equipment

- Connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.
- Consult the manufacturer or field service technician for help.

5 Warranty and Service

Limited Warranty

Instrumentation purchased new from Welch Allyn is warranted to be free from defects in material and workmanship under normal use and service for a period of one year from the date of first shipment from Welch Allyn. This warranty shall be fulfilled by Welch Allyn or its authorized representative repairing or replacing, at Welch Allyn's discretion, any such defect free of charge for parts and labor.

Welch Allyn should be notified via telephone of any defective product and the item should be immediately returned, securely packaged and postage prepaid to Welch Allyn. Loss or damage in shipment shall be at purchasers' risk.

Welch Allyn will not be responsible for loss associated with the use of any Welch Allyn product that (1) has had the serial number defaced, (2) has been repaired by anyone other than an authorized Welch Allyn Service Representative, (3) has been altered, or (4) has been used in any manner other than in accordance with instructions.

THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ANY IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR PARTICULAR PURPOSE, OR OTHER WARRANTY OF QUALITY, WHETHER EXPRESS OR IMPLIED, Welch Allyn WILL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES.

Service

Before returning a product for repair or annual calibration you must obtain authorization from Welch Allyn. An RMA (Return Merchandise Authorization) number is given to you from our service personnel. Be sure to note this number on the outside of your shipping box. Returns without an RMA number will not be accepted for delivery. Please contact your customer service representative at the appropriate Welch Allyn service center listed below for further information.

All equipment used in calibrating the 9600 Plus Calibration Tester is traceable to NIST.

Welch Allyn assumes no responsibility if the 9600 Plus Calibration Tester is calibrated at an unauthorized service center. Send the 9600 Plus Calibration Tester to any authorized service center for recalibration.

To find a service center near you, go to <http://www.welchallyn.com/en/about-us/locations.html>.

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Welch Allyn, Inc.