Welch Allyn® fiber-optic laryngoscope handles

Instructions for use

REF 775243 DIR 80029931 Ver. B Revision date: 2022-10

901087 INSTRUMENT HANDLE

Rx ONLY

Fiberoptik-Laryngoskop-Griffe von Welch Allyn
 Manches pour laryngoscope à fibre optique Welch Allyn
 Welch Allyn-laryngoscoop-handgrepen met glasvezelverlichting
 Impugnatura del laringoscopio a fibra ottica Welch Allyn
 Fiberoptiske laryngoskopohåndtag fra Welch Allyn
 Λαβές λαρυγγοσκοπίων οπτικών ινών Welch Allyn
 Mangos de laringoscopio de fibra óptica de Welch Allyn
 Pegas de fibra ótica para laringoscópio da Welch Allyn
 Welch Allynin kuituoptiset laryngoskopohandrat
 Welch Allyn fiberoptiska laryngoskopohandtag
 Uchwyty światłowodowe do laryngoskopu Welch Allyn
 Rukojeť laryngoskopů s optickým vláknem Welch Allyn
 Welch Allyn fiberoptiske laryngoskopohåndtak
 Welch Allyn® fiberoptik laryngoskop kollar
 Волоконно-оптические рукоятки ларингоскопа Welch Allyn
 Welch Allyn опtotалшықты ларингоскоп тұтқалары


© 2021 Welch Allyn, Inc. All rights reserved. No one is permitted to reproduce or duplicate, in any form, this instructions for use or any part thereof without permission from Welch Allyn. Welch Allyn is a trademark of Welch Allyn, Inc. Hillrom is a trademark of Hill-Rom Services, Inc. Welch Allyn, Inc. ("Welch Allyn") assumes no responsibility for any injury to anyone that may result from (i) failure to properly use the product in accordance with the instructions, cautions, warnings, or statement of intended use published in this manual, or (ii) any illegal or improper use of the product.

Welch Allyn, Inc.
 4341 State Street Road
 Skaneateles Falls, NY 13153 USA
 hillrom.com
 Welch Allyn, Inc. is a subsidiary of Hill-Rom Holdings, Inc.

and EU IMPORTER
 Welch Allyn Limited
 Navan Business Park, Dublin Road
 Navan, Co. Meath C15 A222
 Ireland

Authorized Australian Sponsor
 Welch Allyn Australia Pty Limited.
 1 Baxter Drive
 Old Toongabbie NSW 2146
 Australia

Authorized Representative for Kazakhstan
 TOO Orthodox Pharm
 Uly Dala Avenue 7/4, apt 136, Nur-Sultan 010000,
 Kazakhstan

ENGLISH

Instructions for use

Intended use

The fiber-optic laryngoscope handle is an accessory used with compatible rigid fiber-optic laryngoscope blades which are used to examine and visualize a patient's upper airway and aid placement of a tracheal tube.

About this document

These instructions for use apply to Welch Allyn reusable fiber-optic laryngoscope handles 60813, 60813-LED, 60814, 60814-LED, 60815, and 60815-LED. Welch Allyn reusable fiber-optic laryngoscope handles may be used with Welch Allyn fiber-optic laryngoscope blades MacIntosh 6906X, English MacIntosh 6921X, and Miller 6806X.

Symbol descriptions

For information on the origin of these symbols, see the Welch Allyn symbols glossary: welchallyn.com/symbolsglossary

Reorder number	Product Identifier
Medical device	Lot code
Manufacturer	Humidity limitation
Temperature limit	Consult Instructions for use
Global Trade Item Number	

Separate collection of Electrical and Electronic Equipment. Do not dispose as unsorted municipal waste.

Rx ONLY Prescription only or "For Use by or on the order of a licensed medical professional"

Authorized Representative in the European Community

WARNING The warning statements in this manual identify conditions or practices that could lead to illness, injury, or death. Warning symbols will appear with a grey background in a black and white document.

CAUTION The caution statements in this manual identify conditions or practices that could result in damage to the equipment or other property, or loss of data.

Warnings

- WARNING** Reusable fiber-optic laryngoscope handles must be reprocessed after each use.
- WARNING** The reprocessing procedure and the equipment and materials described must be followed and conducted by persons trained and familiar with medical device reprocessing.
- WARNING** Consult cleaning and disinfecting agent manufacturer instructions for their proper preparation and use.
- WARNING** Repeated reprocessing may degrade elements of the handle. Follow inspection procedures to assure damage has not occurred to the handle.
- WARNING** High-level disinfection and/or sterilization are not achieved by these methods.
- WARNING** Discard any component that shows evidence of damage or deterioration. Contact Hillrom Technical Support to replace components.
- WARNING** Do not modify this equipment. Any modification of this equipment can lead to patient injury. Any modification of this equipment voids the product warranty.
- WARNING** Personnel shall follow their facility policies and procedures and wear appropriate personal protective equipment when handling potentially contaminated equipment.
- WARNING** Laryngoscope equipment is not suitable for use in intense magnetic fields.
- WARNING** Lamps, if left illuminated, could generate sufficient heat to cause burns.

Cautions

- CAUTION** Only use lamp specified. Failure to follow these instructions may cause damage or poor performance of the handle.
- CAUTION** If the device will be unused for several months or longer, remove the batteries prior to storing the device.

Residual risk

This product complies with relevant electro-magnetic interference, mechanical safety, performance, and biocompatibility standards. However, the product cannot completely eliminate potential patient or user harm from the following:

- Harm or device damage associated with electro-magnetic hazards,
- Harm from mechanical hazards,
- Harm from device, function, or parameter unavailability,
- Harm from misuse error, such as inadequate cleaning, and/or
- Harm from device exposure to biological triggers that may result in a severe systemic allergic reaction.

Reprocessing

Reprocessing refers to procedures for cleaning and intermediate-level disinfection. Fiber-optic laryngoscope handles must be reprocessed before first use and between each use according to the following method outlined in this document:

- Cleaning and intermediate-level disinfection

Welch Allyn has validated these instructions as being capable of preparing the handles for re-use. You must ensure that the reprocessing as performed by your facility's personnel, using the facility's equipment and materials, achieves the desired results. This may require validation and routine monitoring of your

facility's actual process. This also may require validation and routine monitoring of the user's actual process. These handles are suitable for 200 reprocessing cycles.

NOTE The main handle and bottom cap components are compatible with the disinfection solution and autoclave methods identified which are provided for facilities who wish to perform either method after cleaning and intermediate-level disinfection.

Cleaning and intermediate-level disinfection

Point of use

1. Separate blade assembly from handle and place handle into suitable containment for subsequent reprocessing. See Figure 1. Do not place handle with sharp devices.
2. Prevent the handle from drying (i.e., wrap/cover in moist germicidal wipe).

Preparation for decontamination

1. Select an appropriate quaternary ammonium isopropanol-based germicidal cleaner labeled suitable for use on healthcare equipment and capable intermediate-level disinfection. Reference EPA-registered disinfectants: www.epa.gov/pesticide-registration/selected-epa-registered-disinfectants. Outside of the U.S., consult applicable regulatory body for equivalent quaternary ammonium isopropanol germicidal cleaner.
2. Remove batteries and lamp cartridge. See Figure 2.

Cleaning and intermediate-level disinfection

1. Follow the germicidal wipe manufacturer's instructions to clean all exposed surfaces of the main handle, end cap, and lamp cartridge.
2. If necessary, brush with a dry, soft-bristled brush and re-wipe to loosen/remove excessive visible soil.
3. After all visible soil is removed, re-wipe to wet all surfaces and allow adequate contact time for disinfection as directed by the germicidal wipe manufacturer.

CAUTION Only use quaternary ammonium isopropanol-based germicidal wipes.

Drying

Allow components to air dry.

Maintenance, inspection, and testing

1. Inspect each component area for damage or deterioration. Inspect (1) handle REF number, (2) main handle, (3) lamp cartridge assembly, (4) batteries, and (5) bottom cap. See Figure 3.

WARNING Discard any component that shows evidence of damage or deterioration. Contact Hillrom Technical Support to replace components.

2. Reassemble handle with new or batteries in known good condition. See Figure 4.

3. Attach handle to a clean and disinfected test blade in known working condition. Verify the following:

- Blade assembly engages and locks onto handle.
- Blade assembly deploys into its locked position on handle AND lamp illuminates.
- Light output is satisfactory.

If the lamp fails to light or output is low, check or replace the lamp and/or batteries.

Storage

Store handle per facility practice to allow device to remain clean, dry, and ready for service.

End of reprocessing instructions for intermediate-level cleaning.

Cold solution disinfection

NOTE The main handle and bottom cap are compatible with this cold solution disinfection method which is provided for facilities who wish to disinfect with solution after cleaning and intermediate-level disinfection.

Disassembly

1. Disassemble the handle and remove batteries and lamp cartridge. See Figure 2.
2. Set batteries and lamp cartridge aside.

Preparation for decontamination

1. Select a 14 day (2.4 – 2.6%) glutaraldehyde disinfectant.
2. Prepare disinfection solution per manufacturer instructions.

Cold solution disinfection

1. Immerse main handle and bottom cap (only) in disinfection solution, for a time duration identified by the disinfection solution manufacturer.
2. Thoroughly rinse all components in potable water, softened water, or per disinfection solution manufacturer instructions to remove disinfecting solution.

Drying

Dry all components with a clean cloth or allow to air dry.

Maintenance, inspection, and testing

1. Inspect each component area for damage or deterioration. See Figure 3.

WARNING Discard any component that shows evidence of damage or deterioration. Contact Hillrom Technical Support to replace components.

2. Reassemble with new or batteries in known good condition. See Figure 4.

3. Attach handle to a clean and disinfected test blade assembly in known working condition. Verify the following:

- Blade assembly engages and locks onto handle.
- Blade assembly deploys into its locked position on handle AND lamp illuminates.
- Light output is satisfactory

If the lamp fails to light or output is low, check or replace the lamp and/or batteries.

Storage

Store handle per facility practice to allow device to remain clean, dry, and ready for service.

Autoclave

NOTE The main handle and bottom cap are compatible with these autoclave methods which are provided for facilities who wish to autoclave after cleaning and intermediate-level disinfection.

Disassembly

1. Disassemble the handle and remove batteries and lamp cartridge. See Figure 2.

2. Set batteries and lamp cartridge aside.

After battery and lamp cartridge removal, select ONE of the following autoclave methods below for the main handle and bottom cap.

Gravity autoclave

Follow equipment manufacturer and facility procedures in the setup and operation of autoclave equipment. Gravity autoclave settings are as follows:

- Temperature: 132 °C (270 °F)
- Exposure time: 3 minutes (unwrapped)
- Minimum dry time: 1 minute

Pre-vacuum autoclave

Follow equipment manufacturer and facility procedures in the setup and operation of autoclave equipment. Pre-vacuum autoclave settings are as follows:

- Temperature: 132 °C (270 °F)
- Exposure time: 3 minutes (unwrapped)
- Minimum dry time: 1 minute

Maintenance, inspection, and testing

1. Inspect each component area for damage or deterioration. See Figure 3.

WARNING Discard any component that shows evidence of damage or deterioration. Contact Hillrom Technical Support to replace components.

2. Reassemble with new or batteries in known good condition. See Figure 4.

3. Attach handle to a clean and disinfected test blade assembly in known working condition. Verify the following:

- Blade assembly engages and locks onto handle.
- Blade assembly deploys into its locked position on handle AND lamp illuminates.
- Light output is satisfactory.

If the lamp fails to light or output is low, check or replace the lamp and/or batteries.

Storage

Store handle per facility practice to allow device to remain clean, dry, and ready for service.

Maintenance

Replace the lamp

WARNING Lamps, if left illuminated, could generate sufficient heat to cause burns.

CAUTION Use only Welch Allyn replacement lamps

REF 06000-LED to ensure proper illumination alignment.

1. Unscrew bottom cap of handle clockwise and remove batteries. See Figure 2.
2. Remove lamp cartridge assembly from main handle by applying finger pressure in the direction shown by the arrow. See Figure 2.
3. Remove shroud cap from shroud by rotating shroud cap counterclockwise (1) shroud cap (2) lamp, (3) shroud. See Figures 5 and 6.
4. Rotate lamp cap counterclockwise to remove. See Figure 6.
5. Replace lamp and rotate shroud cap clockwise to tighten. See Figure 7.
6. To replace lamp cartridge assembly in the main handle, invert handle, then gently slide the cartridge down the inside of handle, tipping it side to side until the shroud cap exits opening on top.
7. Re-insert batteries and apply slight pressure to set the cartridge in place.
8. Replace and tighten bottom cap.
9. Reprocess repaired assembly as appropriate per these instructions.

Replace the batteries

1. Unscrew bottom cap of handle and remove batteries. See Figure 2.
2. Insert batteries and reinstall bottom cap. See Figure 4.

NOTE Alkaline batteries are supplied with your handle for maximum performance and are recommended as replacements; however carbon-zinc batteries may also be used.

- Medium handle, 60813 model, use two "C" size
- Penlight handle, 60814 model, use two "AA" size
- Stubby handle, 60815 model, use two "AA" size

3. Verify lamp and blade engagement/operation using a known working test blade.
4. Reprocess repaired assembly as appropriate per these instructions.

Disposal

Users must adhere to all federal, state, regional, and/or local laws and regulations as they pertain to the safe disposal of medical devices and accessories. If in doubt, the user of the device should first contact Hillrom Technical Support for guidance on safe disposal protocols.

PATENT/PATENTS

hillrom.com/patents

May be covered by one or more patents. See above Internet address. The Hill-Rom companies are the proprietors of European, US, and other patents and pending patent applications.

Hillrom Technical Support

For information about any Hillrom product, contact Hillrom Technical Support at hillrom.com/en-us/about-us/locations.

Notice to users and/or patients in the EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Specifications

Temperature limits

Operating: 32° F to 104° F (0° C to 40° C)
 Storage/Transport: -4° F to 120° F (-20° C to 49° C)

Standards and compliance

The device complies with the following standards: ISO 7376, IEC 60601-1, IEC 60601-1-2

Country-specific standards are included in the applicable Declaration of Conformity.

Lot code

For lot code YY-JJJ, YY= Last 2 digits of the year, and JJJ=consecutive day of the Julian year.

Guidance and manufacturer's declaration

For information about electromagnetic compatibility (EMC), see the Hillrom website:

welchallyn.com/en/service-support/emc.html

A printed copy of EMC Emissions and immunity information can be ordered from Hillrom for delivery within 7 calendar days.

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