



WELCH ALLYN MULTI-PARAMETER VITAL SIGNS DEVICES - SCANNER OFFERINGS



SCANNER FUNCTIONALITY	CONNEX SPOT MONITOR	CONNEX VITAL SIGNS MONITOR	CONNEX INTEGRATED WALL SYSTEM
Barcode Only	7000-916HS* 7000-CS6080***	6000-916HS**	6000-916HS**
High Frequency RFID + Barcode	7000-916HSR*	6000-916HSR**	6000-916HSR**
Low Frequency RFID + High Frequency RFID + Barcode	7000-916HS1RS*	7000-916HS1RS**	7000-916HS1RS**

COIL CORD REPLACEMENT KIT COMPATIBILITY	6000-916HS	7000-916HS	7000-CS6080	6000-916HSR	7000-916HSR	7000-916HS1RS
JADAK Scanner Coiled Cord Replacement Kit (106594)	✓	✓	✗	✓	✓	✓
Zebra Scanner Coiled Cord Replacement Kit (108782)	✗	✗	✓	✗	✗	✗

SCANNER HOLSTER REPLACEMENT KIT COMPATIBILITY	6000-916HS	7000-916HS	7000-CS6080	6000-916HSR	7000-916HSR	7000-916HS1RS
JADAK Scanner Holster Kit (7000-HOLSTER)	✓	✓	✗	✓	✓	✓
Universal Scanner Holster Kit (108781)	✓	✓	✓	✓	✓	✓

*Compatible with **Connex** Spot Monitors (CSM) with 1.30.00 firmware or later.

Compatible with **Connex Vital Signs Monitors (CVSM) and **Connex** Integrated Wall Systems (CIWS) with 2.42 firmware or later.

***Compatible with **Connex** Spot Monitors (CSM) with 1.40.10 A0002 firmware or later.



As the 7000-916HS, 7000-916HSR, 6000-916HS, 6000-916HSR and 7000-916HS1RS scanners share similar housings which can make it difficult to differentiate them, you can distinguish between them using their model numbers which can be found at the bottom on the back of the scanner housing.

For more information, please contact your Baxter representative.

Baxter.com | Baxter International Inc.

Rx only. For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

Baxter, Connex and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.

JADAK is a trademark of Novanta Corporation. Zebra is a trademark of Zebra Technologies Corp.

US-FLC45-240020 (v1.0) 06/2024